

**Dr Anne Sophie MARTY**

Thiers Ophthalmology Practice, Thiers Ophthalmology Clinic,  
Bordeaux, France

**Key words**

Dry Eye Syndrome (DES) - Dry Eye - Ocular Surface Disease (OSD) -  
Meibomian Gland Dysfunction (MGD) - Tear film- Refractive surgery -  
eTAO – C.DIAG® - Meibography - C.STIM® - Intense Pulsed Light (IPL)

**Summary**

**Diagnosis and treatment** of a female patient with mixed **Dry Eye Syndrome (DES)** aggravated by **refractive surgery**.

A full diagnosis was made via clinical examination and examinations using the **C.Diag®** imaging platform.

Symptomatic treatment combined with **Intense Pulsed Light** treatment using the **C.Stim®** IPL system was started for the patient. **After 3 C.Stim® IPL sessions, the patient's quality of life was restored** due to a clear improvement in symptoms and an improvement in MGD, with better meibum expression and quality.

**1**

Clinical  
examination

**2**

Para-clinical  
examination

**3**


C.Stim® IPL  
treatment

**4**

Results  
at 2 months

**Anamnesis**

- 30-year-old female patient
- Deputy director in the medico-social sector
- Laser eye surgery (LASIK) in 2018
- Symptoms since 2019
- Dry eye sensation, stinging, watering, burning, redness, nighttime pain (OD>OS).
- **Symptoms frequency: 5/10**
- **Symptoms intensity: 8/10**
- **DEQ-5=13**



# 1

## Clinical examination

### Existing treatment:

- Eyelid care, blinking exercises
- Artificial tears as required (Hyloconfort®, Hylovis lipo®, Liposic® or Vitamin A® ointment, Naabak®)
- Intolerance to Azyter®, Ikervis®, Doxycycline®
- Lacrimal plug in left lower punctum

### DES Risk factors:


- LASIK surgery
- Screen use
- Cutaneous rosacea
- History of chalazions
- Previous contact lens use
- Dust mite and pollen allergy and long-term antihistamine use.

### Ophthalmologic exams:

OD		OS
10/10 P2	AV	10/10 P2
+0,25 (-0,25 110°)	Refraction	plan
12	IOP	15

### Slit lamp exams:

- No reveal superficial punctate keratitis
- Fluorescein-filling defect (presence of spots)
- BUT average of 7 sec//8 sec
- Lid margin telangiectasias
- Thick meibum
- Conjunctival hyperemia
- Conjunctival papillae
- LIPCOF




# 2

## Para-clinical examination


### 1. Epithelial mapping:

- Onset of concentric epithelial atrophy



### 2. eTAO score

- eTAO exam report ables to measure the DES severity score



### 3. Ocular surface analysis with C.DIAG®


#### Tear meniscus height:

lower OD= 0.17 mm // OS = 0.23 mm

#### Interferometry: correct ODS

**NIBUT:** reduced to 8 sec OD / 10.4 sec OS


→ **Slight tear instability**



→ **Aqueous deficiency**

#### Transillumination meibography:

stage 1 MG atrophy (MG loss < 25%),  
stage 2 telangiectasias, note the presence  
of the OS lacrimal plug (arrow)



**Blink analysis:** pathological OD>OS

- **Insufficient blink rate in the right eye**
- **More than 50% incomplete blinks**



## Focus on contact lens intolerance

=> A patient who used to wear contact lenses but became intolerant to them over time should raise concerns of developing MGD.

=> Contact lens intolerance is a frequent cause for consultation with a view to refractive surgery.

=> Wearing contact lenses: - can lead to abnormal blinking,



- splits the tear film in two (pre- and post-lens),
- makes the lipid layer thinner and more uneven (lipid layer),
- causes corneal hypoesthesia leading to a lower tear turnover rate (aqueous layer)
- reduces the number of mucus cells (mucin layer).

=> Screen work is the main cause of incomplete blinking.

**Wearing contact lenses while working on a screen is often uncomfortable due to the resulting MGD.**


## Diagnosis

Mixed dry eye due to moderate Meibomian gland dysfunction (MGD), along with post-refractive surgery aqueous deficiency

→ **An ophthalmological assessment was carried out to identify systemic dry eye syndrome prior to refractive surgery.**

## Initial treatment

- Lacrimal plugs - OD
- Eyelid care and blinking exercises
- Systane Ultra® mucous membrane protector
- Suggestion of an allergen immunotherapy by an allergy specialist



## 3 month check-up

Symptoms are reduced but are still present despite the patient complying with the treatment.

- ▶ **C.STIM® IPL treatment is initiated**
- ▶ **Continuation of symptomatic treatment**

### 3 C.STIM® IPL treatment

- Treatment protocol: 3 sessions on D0, D15 and D30
- 4 shots per side per session at a fluence of 8 J/cm<sup>2</sup>
- Protective goggles for both the patient and the doctor
- Meibum expression using forceps done after each session



### 4 Results at 2 months

- **Clear improvement in symptoms experienced:**
    - ▶ **Feedback from the last session:** « eyes are more comfortable ».
    - ▶ **Artificial tears:** she keeps putting them on the right only during the day and uses Liposic® before bed.
    - ▶ **Warming mask and massaging:** she performs eyelid hygiene 3 times a week.
    - ▶ **Blinking exercise:** she performs them during screen breaks at work.
    - ▶ **Environment:** she avoids heating the bedroom.
  - **Clinical examination:**
    - ▶ No SPK
    - ▶ BUT 10 sec//10 sec
    - ▶ Decrease of lid margin telangiectasias
    - ▶ More fluid meibum OS>OD during meibum expression
- ➔ **Symptomatic treatment must be continued over the long term**
- ➔ **IPL treatment repetition if necessary**

## Conclusions

- C.STIM® IPL treatment is **rapid, effective and safe**.
- **Dry eye** following **refractive corneal (or crystalline lens) surgery** is a **major cause of post-surgical dissatisfaction**.
- **IPL treatment** can easily be offered to patients already suffering from **pre-operative and/or post-operative MGD** when local symptomatic treatment is not effective enough.