



Dr Marie-Caroline TRONE

Hospital-based Ophthalmologist
Saint-Étienne University Hospital – France

Key words

Dry eye syndrome (DES) – Meibomian gland dysfunction (MGD)
Dry eye – Blepharitis – Demodex - Ocular surface – Tear film
OSDI – Intense pulsed light (IPL) - C.Stim® – TFOS-DEW II

Summary

Diagnosis and treatment of **Dry Eye Syndrome (DES)** in a male patient with **Demodex blepharitis leading to MGD**. A full diagnosis was performed by means of clinical examination and in vivo confocal microscopy. **C.STIM® intense pulsed light treatment** plus eyelid hygiene was initiated in this patient. After three months, **a significant improvement in functional signs and meibum quality was observed, with almost no Demodex sleeves visible**.

1

Clinical examination

2

Confocal microscopy examination

3

IPL treatment C.Stim®

4

Results at three months

Patient history

- 57-year-old man
- Teacher
- Eyes stuck together in the morning upon waking
- Symptoms mainly occur in the morning, with fluctuating blurred vision upon blinking
- Dry eye syndrome that has been debilitating for several months, with the prescribed symptomatic treatment offering little relief: stinging, itching and photophobia



1

Clinical examination

Interview :


- Identification of DES risk factors:
=> History of OD corneal transplant
- OSDI = 35.4
- Good compliance with symptomatic treatment:
artificial tears, eyelid care (massage)
- Free margins cleaned in the evening with tea tree oil wipes: not tolerated
(inflammatory reaction of the eyelids leading to discontinuation after a few days)

Eye examination:

OD	VA	OS
8/10 Pa 2 -2.00 (-2.25 at 110°) add+2.00		10/10 Pa2 -1.50 (-0.75 at 180°) add+2.00
15 mmHg	IOP	17 mmHg
Transplant clear A few spots of SPK lower down	Cornea	Cornea clear A few spots of SPK lower down

Slit lamp examination: Mild to moderate MGD/Demodex blepharitis

- Incomplete blinking
- OD BUT 5 seconds and OS BUT 4 seconds with significant tear instability




Moderate MGD with sleeves around eyelashes and crusting+++

2

In vivo confocal microscopy (Vivascope 3000)

Demodex blepharitis visible under in vivo confocal microscopy: reflectance at an infundibulum of a eyelash filled with elongated, parallel bodies corresponding to *Demodex folliculorum*



Diagnosis

Evaporative dry eye syndrome with meibomian gland dysfunction and Demodex blepharitis


3

C.STIM® I.P.L. Treatment

- Symptomatic treatment continued


- C.STIM® IPL treatment

- ▶ Three sessions on D0, D15 and D45
- ▶ Four shots per side per session at a fluence of 8 J/cm²
- ▶ Protective goggles worn by patient and doctor
- ▶ Meibum expression after each session (with forceps)



4

Results at three months




- **Improvement in functional signs**
- **Improvement in meibum quality** and almost total disappearance of sleeves



Conclusions

IPL treatment is effective against Demodex blepharitis in combination with, or as an alternative to, symptomatic treatment (e.g. cleaning the free margins with tea tree oil wipes).

Focus on Demodex



- Responsible for symptomatic and non-symptomatic blepharitis/MGD
- Sleeves around the roots of the eyelashes
- Sometimes visible under in vivo confocal microscopy
- Conclusive diagnosis: parasitological examination of eyelashes



Figures: SFO [French Society of Ophthalmology] Report 2015, Ocular Surface

Bibliography

1. Parasympathetic Innervation of the Meibomian Glands in Rats - Mark S. LeDoux et al. - *Investigative Ophthalmology & Visual Science*, October 2001, Vol. 42, No. 11
2. Characterization of the innervation of the meibomian glands in humans, rats and mice - Bründl, M. et al. *Annals of Anatomy* (2021), Vol. 233.
3. Neurotransmitter Influence on Human Meibomian Gland Epithelial Cells - Wendy R. Kam and David A. Sullivan - *Investigative Ophthalmology & Visual Science*, November 2011, Vol. 52, No. 12
4. The Dopaminergic Neuronal System Regulates the Inflammatory Status of Mouse Lacrimal Glands in Dry Eye Disease - Ji, Yong Woo et al. *Investigative Ophthalmology & Visual Science* (2021), Vol. 62.
5. TFOS DEWS II pathophysiology report - Anthony J. Bron, et al. - *The Ocular Surface*, 2017, p 441 to 515
6. The neurobiology of the meibomian glands - Cox SM, Nichols JJ - *Ocular Surface*, July 2014
7. Multicenter Study of Intense Pulsed Light Therapy for Patients with Refractory Meibomian Gland Dysfunction - Reiko Arita, et al. - *Cornea* Volume 37, Number 12, December 2018
8. Rosacea: Molecular Mechanisms and Management of a Chronic Cutaneous Inflammatory Condition - Yu Ri Woo, et al. *International Journal of Molecular Sciences*, September 2016
9. Rosacea: Epidemiology, pathogenesis, and treatment - Barbara M. Rainer et al. - *DERMATO-ENDOCRINOLOGY* 2018, VOL. 9, NO. 1, e1361574 (10 pages)
10. Treatment of ocular rosacea - Edward Wladis et al. - *Survey of Ophthalmology* (2018), Vol.63.
11. Improved telangiectasia and reduced recurrence rate of rosacea after treatment with 540 nm-wavelength intense pulsed light: A prospective randomized controlled trial with a 2-year follow-up - Luo, Y. et al. - *Experimental and Therapeutic Medicine* (2020), Vol. 19.
12. Therapeutic Effect of Intense Pulsed Light on Ocular Demodicosis - Zhang, X., et al. - *Current Eye Research* 2019, Vol. 3.
13. Intense Pulsed Light Therapy for Patients with Meibomian Gland Dysfunction and Ocular Demodex Infestation - Cheng et al. - *Current Medical Sciences* (2019), Vol.39.
14. Long-term effects of intense pulsed light treatment on the ocular surface in patients with rosacea-associated meibomian gland dysfunction - Seo Kyoung Yul et al. - *Contact Lens and Anterior Eye* (2018), Vol. 41.
15. TFOS DEWS II Tear Film Report - Willcox Mark et al. - *The Ocular Surface* (2017), Vol.15.
16. Intense Pulsed Light for the Treatment of Dry Eye Owing to Meibomian Gland Dysfunction - Vigo, L. et al. - *Journal of Visualized Experiment* (2019), N°146.
17. Meibum Expressibility Improvement as a Therapeutic Target of Intense Pulsed Light Treatment in Meibomian Gland Dysfunction and Its Association with Tear Inflammatory Cytokines - Choi, M. et al. - *Scientific Reports* (2019), Vol.9.
18. TFOS DEWS II pain and sensation report - Belmonte Carlos, et al. - *The Ocular Surface* (2017), Vol.15.
19. Analysis of Cytokine Levels in Tears and Clinical Correlations After Intense Pulsed Light Treating Meibomian Gland Dysfunction - LIU, R et al. - *American Journal of Ophthalmology* (2017).
20. Effect of inflammation on lacrimal gland function - Driss Zoukhri - *Experimental Eye Research*, May 2006; 82(5): 885-898