

Compact Touch[®]



**COMPACT IN
DESIGN,
DEFINITIVE IN
DIAGNOSIS**




A/B/P Ultrasound Platform

Compact Touch[®]


■ QUANTEL MEDICAL REDEFINES THE STANDARD IN B-SCAN IMAGING

The Compact Touch[®] benefits from a **new generation of 15 MHz B probe** with a **resolution increased by 30%**. It allows for a better visualization of the eye structures and the orbit hence a better diagnosis.


Small in size, this probe benefits from an **excellent ergonomics to facilitate its handling and use**.




**LUXATED
LENS**



**RETINAL
DETACHMENT**




**MALIGNANT
MELANOMA**



**SUBTLE AV
FISTULA**

■ BIOMETRY IN B MODE

Compact Touch® has the **exclusive technology : Biometry in B-Mode**, it allows the **measurement of axial length**, just in the macula, especially in eyes of unusual geometry, like patients with long myopic eyes associated with staphylomas, while ruling out any other pathology of the posterior pole.



HIGH MYOPIA

■ BIOMETRY AND IOL CALCULATION

Ultrasound biometry is the **only technology** suitable for axial measurement of all eye types, regardless of cataract density and media opacity. Ultrasound axial length measurement offers similar levels of precision to optical measurement (0.03 mm with immersion technique)¹.


The IOL calculation function can compare the results of different types of IOLs and formulas. A total of **12 calculation formulas** are available, including formulas for post-refractive surgery patients. IOLs can be calculated in intervals up to **0.25D**.

■ PROBEAM™ BIOMETRY PROBE*

ProBeam™ is exclusive to Quantel Medical. The probe generates a laser beam that creates a focal point for the patient to look at with the other eye macula: this facilitates measurement while increasing the precision of the probe².




Quantel Medical, the world leader in ophthalmic ultrasound



■ A SLEEK DESIGN, WITH INCREASED ERGONOMICS

A weight **less than 4 kg**, Compact Touch® is aptly named, since its dimensions have been reduced to suit for all types of clinical settings.

- It is more easily transportable, thanks to its foldable and **reclining carrying handle**.
- Fanless, it is as **silent** as a tablet.
- Its **VESA fixation system** offers the possibility to be **mounted either on a wall, on an articulated arm or on a mobile cart**, to enhance ergonomics around the patient*.



■ PACHYMETRY*

Pachymetry is essential for **refractive surgery and glaucoma diagnosis**. The Compact Touch® has several measurement modes, **offering a precision of ± 5 microns** and a **measurement range of 200 to 999 microns**.

IOP measurements can be corrected using integrated tables of correlations between IOP and corneal thickness, including Ehlers, Doughty and Dresden.

(*) Option



■ A CONNECTED ULTRASOUND PLATFORM

- With **DICOM interface**, the Compact Touch® can import (worklist function) and export (storage function) images and patients reports to the PACS. Reports and images printout is also possible either on a **DICOM printer**, or a **local printer**.
- **Video sequences** (Cineloops) can be sent in **DICOM format**.
- For easier use, a **wireless keyboard** and **mouse** can be connected.



B MODE

Grey levels:	256
Adjustable gain:	20 to 110 dB
Time Gain Control (TGC):	0 to 30 dB
Dynamic range adjustment:	25 to 90 dB
Storage of still images and video sequences (up to 40 second duration) on hard disk	
Image post-processing tools:	calipers, areas, markers, comment

15 MHz probe

Transducer frequency:	15 MHz
Angle of exploration:	50°
Depth of exploration:	60 mm (2.36")
Focus:	24 mm (0.94")
Axial resolution:	115 µm
Lateral resolution:	400 µm
Frame Rate Acquisition:	up to 16 Hz

BIOMETRY

Adjustable gain:	20 to 110 dB
Time Gain Control (TGC):	0 to 30 dB

11 MHz Probe

Transducer frequency:	11 MHz
Tip diameter:	7 mm (0.28")
Electronic resolution:	0.03 mm (0.002")
Depth:	60 mm (2.4") on 1536 points
Contact and immersion techniques compatible	
Aiming beam:	LED or laser pointer ProBeam™

Axial length measurements

Ultrasound propagation velocity adjustable per segment (anterior chamber, lens, vitreous) and IOL and vitreous material
 Built-in pattern recognition: phakic, aphakic, PMMA, acrylic and silicone material for pseudo-phakic eye types
 Automatic calculation of standard deviation and average total length (series of 10 measurements)
 Acquisition modes: automatic, auto + save, manual
 Automatic detection of scleral spike

IOL calculation

SRK-T, SRK-II, HOLLADAY, BINKHORST-II, HOFFER-Q, HAIGIS
 Post-op refractive calculation:
 - Pre-op and Post-op refraction, Pre-op and Post-op keratometry
 - 6 different methods for keratometric correction and implant calculation:
 History derived, refraction derived, contact lens method, Rosa regression, Shammas regression, Double K/SRK-T (Dr. Aramberti's formula)
 9 values bracketed for desired ametropia for each IOL (IOL increment steps: 0.25D or 0.50D). Simultaneous display of 4 different IOL calculations

DATA MANAGEMENT

Built-in physician and patient database
 Exportation of still images and video sequences
 Customisable digital and printed report
 DICOM compatible (Worklist, Storage, Print)*
 EMR compatible
 Compatible with PC and USB video printers

TECHNICAL SPECIFICATIONS

PACHYMETRY*

Transducer frequency:	20 MHz
Tip diameter:	1.2 mm (0.05")
Method:	contact
Convergence:	0.5 mm (0.02") from the tip
Angle:	45°

Corneal thickness measurements

Measurement range:	200 to 999 microns
Number of measurements:	1 to 10
Precision:	± 5 microns
Velocity:	adjustable
Methods:	central measurement or cartographic map (automatic, continuous, scanning)
Cartographic map:	user - 9C8L - 9C4L - 5C8L - 5C4L - 9C - 5C - 8L - 4L

I.O.P. correlation tables

Tables correlating intraocular eye pressure and corneal thickness:
 Ehlers + Doughty + Dresdner + unlimited user-defined tables

Specifications

Bias correction:	up to 120%
------------------	------------

GENERAL INFORMATION

Back-lit LCD colour touch screen monitor (resolution 1024 x 768 px)

Electrical requirements

Power supply:	100-240 Vac ±10% single phase with earth
Frequency:	50/60 Hz
Power:	60 W max

Features

Overall dimensions:	26.8 cm (W) x 4.0 cm (D) x 24.6 cm (H) 10.6" x (W) x 1.6" (D) x 9.7" (H)
Touch screen dimensions:	21 cm (W) x 16 cm (H) - 8.3" (W) x 6.3" (H)
Weight:	3.5 kg (7.7 lbs)
Ports:	4 USB, 1 ethernet, 1 HDMI

Peripherals and accessories included in the basic configuration

Footswitch
 Bluetooth mouse

Peripherals and accessories in option*

External PC printer Windows Operating System compatible (USB or Wifi)
 Video printer with USB connection

(* Option

Specifications are subject to change without notice.
 Non contractual pictures

BIBLIOGRAPHY

1. Comparison of immersion ultrasound biometry and partial coherence interferometry for intraocular lens calculation according to Haigis – W. Haigis et al. – Graefes Arch Clin Exp Ophthalmol. 2000 Sep
2. New laser fixation device for ultrasound biometry – M. Charles – Ophthalmol. Clin. Exp. 2007



www.lumibird-medical.com

Manufacturer

Quantel Medical
 1 Rue du Bois Joli – CS40015
 63808 Cournon d'Auvergne – FRANCE
 Tel.: +33 (0)4 73 745 745
 Email: contact@quantelmedical.fr
 ISO 9001 : 2015 – ISO 13485 : 2016

Headquarters

Lumibird Medical
 1 Rue du Bois Joli – CS40015
 63808 Cournon d'Auvergne – FRANCE
 Tel.: +33 (0)4 73 745 745



QUANTEL MEDICAL - ELLEX - OPTOTEK MEDICAL

© 2022. Quantel Medical, Compact Touch[®] and ProBeam™ are trademarks of Quantel Medical. All rights reserved.