

# DRY EYE

## QUANTEL MEDICAL

## DIAGNOSIS

NEW


**c.diag**

### DIAGNOSTIC PLATFORM

#### DRY EYE / OCULAR SURFACE DISEASES



- **Fast & full automatic exams with one click complete report:** dry eye qualification, quantification and classification  
*Blinks, tear stability, lipid layer, tear meniscus, meibography, eyelids and ocular surface (blepharitis, Demodex, bulbar redness, etc.), exams in fluorescein*
- **Outstanding image quality,** standard infrared for quantification and exclusive transillumination infrared for qualification
- **Robust AI algorithms:** high precision and reliable results
- **Intuitive handling:** integrated design and friendly software



Not FDA Approved



**LEARN MORE**

## TREATMENT

**c.stim**

### IPL SYSTEM, INTENSE PULSED LIGHT



- **Treatment of the root causes of dry eye:** MGD, inflammation, Demodex and neuropathic pain
- **Quick and easy treatment (<10 min)**
- **Safe and efficient treatment** thanks to the unique handpiece with an optical tip incorporating the **stim-ULI™ technology** (Uniform Light Intensity) and the train of **5 regulated pulses**.
- **Treatment on the majority of skin types** thanks to the light spectrum, the train of 5 regulated pulses and to the water cooling of the handpiece...

Not FDA Approved



**LEARN MORE**



**LUMIBIRD**  
MEDICAL

# DRY EYE

COMPLETE OFFER

 **c.suite™**



Find out more about the products  
in the **DRY EYE** range



**LUMIBIRD®**  
MEDICAL