



THE C-EYE DEVICE

BY EMAGINE

CE marked for the following indications: keratoconus, ectasia after LASIK/PRK, pellucid marginal degeneration, infectious keratitis, sterile corneal melting, bullous keratopathy.



Slit Lamp Mount

For modern CXL: fits on Haag Streit and Zeiss design-slit lamps.



Range of Intensities

Delivers 3, 9, 15, 18 and 30 mW/cm².



Table Mount

For classic CXL: table mount accessory available.



Battery-Operated

Charge via USB-C for maximal mobility.



Thickness-Adjusted Profile

Beam profile delivers more energy in the corneal periphery.



Pulsed and Continuous Light

For various keratoconus & infectious keratitis settings.



Preset Protocols

Covering the most widely used CXL settings.

www.emagine-eye.com

The C-Eye™ device and Ribo-Ker™ solution are CE marked.

The C-Eye™ device and Ribo-Ker™ solution are not available for sale in the United States.

ONE DEVICE TO PERFORM BOTH...



CXL at the slit lamp
The modern way



CXL in the operating room
The classic way - with a table mount



Innovative

The C-Eye™ device brings cross-linking technology to the slit lamp. Simple, efficient and safe from the novice to the expert user.



Versatile

Perform CXL at the slit lamp. Or use the C-Eye device just like any other CXL device: on a table mount, in your operating room.



Portable

Work in different places? Take the C-Eye™ device wherever you go. Charge it via USB-C, or even use a powerbank on-the-go.

RIBO-KER™

ONE RIBOFLAVIN FOR KERATOCONUS AND KERATITIS TREATMENT

EMAGine is proud to offer Ribo-Ker™, its exclusive CE-marked riboflavin.

No more confusion with different types of riboflavin: EMAGINE has developed one solution for all applications.

- ✓ 0.1% hypoosmolaric riboflavin
- ✓ CE marked
- ✓ Ready for all epi-off CXL and all PACK-CXL applications
- ✓ Can be used on normal thickness and thin corneas
- ✓ Comes in a procedure kit with a speculum and a consumable cap for the C-Eye device



www.emagine-eye.com

The C-Eye™ device and Ribo-Ker™ solution are CE marked.
The C-Eye™ device and Ribo-Ker™ solution are not available for sale in the United States.